



STEP IN BATH & SHOWER  
CONVERSION SYSTEM

## TAKING CARE OF YOUR CAP

### Storing your cap

Your cap can be left in place when not in use or it can be stood up on its end in a convenient location out of the way.

It is recommended that you leave the handles in the open position if you are leaving the cap on the bathtub when not in use.

### Cleaning your cap

The cap can be cleaned with warm soapy water. Use a mild liquid detergent such as dish washing detergent to clean the cap. Rinse with clean water. For stubborn stains or scuff marks, we suggest using a Mr. Clean® Magic Eraser®.

DO NOT use any chemicals; they could adversely affect the finish or function of the cap, thereby voiding your warranty.

The seal is made of silicone. All you should require is a wet cloth to clean off any visible soap scum or dirt. Do not use sharp instruments or pull on the seal as there is a risk of damaging the seal.

Handle assemblies are made of chrome plated castings. Warm water should be sufficient for cleaning them. They can be buffed with a soft

cloth to bring out the shine. It is not recommended that any form of chrome cleaner be used as it may transfer to the main panels and adversely affect them.

The remaining hardware is stainless steel or zinc plated steel. All hardware should be dried following use to ensure that it does not rust.

### Yearly check of the cap and walk through adhesive seal

It is strongly recommended that you check the adhesive seal between the walk through insert and bathtub as well as water testing your cap on an annual basis.

Any unusual water loss when bathing could signify that there is a leak present. In this case you should conduct a water test or contact your installer.

**Please be sure to read the accompanying warranty information in order to understand what is covered and what your responsibilities are in maintaining your warranty coverage. Please note that it is the homeowner's responsibility to ensure that the insert adhesive seal and the cap seal integrity is maintained.**

### IMPORTANT SAFETY INSTRUCTIONS

- ❑ **WARNING: Risk of personal injury.**  
To avoid injury, exercise caution when entering or exiting your walk in tub.
- ❑ **WARNING: Risk of personal injury.**  
To avoid injury, do not use the walk through insert as a "step". Although the cap may support your weight it has not been designed with this purpose in mind.
- ❑ **WARNING: Risk of personal injury.**  
To avoid injury, do not sit or place your weight on the cap. It has not been designed for this. Undue weight could cause the cap to break or crack.
- ❑ **WARNING: Risk of personal injury.**  
Water temperatures over 100° F (38° C) may be hazardous to your health. Check and adjust water temperature for your personal comfort.
- ❑ **WARNING: Risk of personal injury.**  
Risk of electric shock; do not permit the use of electric appliances (such as hair dryer, lamp, telephone, radio, television, etc) within four feet of your bathtub.
- ❑ **WARNING: No food or Alcoholic beverages.**  
Use of your bathtub immediately after meals is not recommended. Avoid alcohol consumption before or during the bathing. Alcoholic beverages can cause drowsiness or hypothermia resulting in loss of consciousness or even drowning.

## SAVE THESE INSTRUCTIONS!

### OBTAINING ADDITIONAL PRODUCT INFORMATION

BATHE AGAIN INC.

800-600-7842 | 416-709-7072 | info@quicktub.com | www.quicktub.com